

2014 Equipment Catalog

# Unchallenged for Quality



**PHSChiropractic**  
by Pivotal Health Solutions


## NEW 2014 Moving Innovation Forward



### TradeFlex Series Built To Be The Industry's Best.

Starting at \$2,815

The TradeFlex table allows the Dr. to perform optimum correction of spinal subluxations or disc problems, yet can still function as an exam, adjusting, or therapy table. This versatility will save the Doctor time and money. This workhorse table is also available with an electric flexion option that offers variable flexion speed and treatment timer.



### ERGOSTYLE FX Long Axis Distraction With Manual or Auto Flexion

Starting at \$7,158 - Page 6

Long axis distraction – administered via hand controller, this feature allows the doctor to administer long axis distracting of the pelvic section up to 5".



### HTG 100P Three Drop Adjustment Table

Starting at \$1,425 - Page 9

Your new HTG 100P table has been designed and engineered to comply with the rigid standards of practitioners and is designed to provide you with years of service. The HTG 100P comes with many standard features, five heights and a 450 lb weight capacity.



### Apollo 2,000 mW Cold Laser A Small Powerful Hand Held Laser

Starting at \$3,999 - Page 14-15

Small, lightweight unit fits in the palm of your hand. Simple on and off button, no programming necessary. Includes: unit docking station, carry case, 2 safety goggles and everything you and your patient need for a successful treatment.



# The All New TradeFlex Series

## Built To Be The Industry's Best

### TradeFlex - E9017 Elevating Manual Flexion Table

Starting at \$2,815

### TradeFlex - E9018 Elevating Auto Flexion Table

Starting at \$5,015

The TradeFlex Series features electric tension controls, optional cervical, thoracic, lumbar, and pelvic drops. Also the ability to perform figure 8 range of motion. The pelvic drop can be used in any position. Introducing optional Thoracic and Lumbar drops in a single chest cushion. P-A and lateral flexion can be locked in any position. The table also features a smooth range of +/- 15 degrees of lateral flexion and a flexion depth of 15 degrees to accommodate almost any angle of drive. A removable flexion control t-bar allows the doctor to move freely around the table which improves patient and doctor comfort. The TradeFlex table allows the Dr. to perform optimum correction of spinal subluxations or disc problems, yet can still function as an exam, adjusting, or therapy table. This versatility will save the doctor time and money. This workhorse table is also available with an electric flexion option that offers variable flexion speed and treatment timer.



E9017 Manual Flexion  
shown with options

[See our website for options available](#)

# The All New TradeMark

## Re-Engineered For Excellence

The new TradeMark table can't be ignored, featuring new aesthetics and NOW available with all THREE DROPS; Thoracic and Lumbar in a single chest cushion, and Pelvic!

### TradeMark - E9012

Starting at \$1,899

The Trademark Stationary is an affordable alternative for a growing practice. With its impressive list of standard features, this table is an excellent choice for a starter table. The Trademark table is a low priced table that offers the same high quality construction as other Pivotal Health Solutions tables.



E9012 TradeMark  
shown with options

[See our website for options available](#)

### TradeMark Features

- Adjustable Face Cushions
- Extension / Flexion Headpiece
- Paper Roller, Cutter & Guide Bar
- 7.5" Ankle Extension
- Break-Away & Elevating Abdominal

# The All New ErgoStyle™ ES2000 Gen II

## World-class Patented DROPS, Upgradeable In The Field

### ErgoStyle ES2000 - ES9300 Elevating Adjusting Treatment Table

Starting at \$4,034

### ErgoStyle ES2000 - ES9320 Elevating Adjusting Treatment Table with 3 Auto-Cocking Drops

Starting at \$6,504

ErgoStyle™ ES2000 changes how all adjusting tables are viewed. The ES2000 is ergonomically designed to meet the physical needs of the clinician while providing patients with the highest level of comfort as they move from position to position. The modular design of the ES2000 allows for quick modifications and upgrades in the field. Upgrades can be added later as your techniques and practice grow.

### ErgoStyle ES2000 Features

- Electric Elevation from 20" to 31" with single foot control and expanding covers over lifting mechanism
- Tilting headpiece with adjustable face cushions
- Stationary Chest
- Fixed Pelvic Section
- Dual Auto-Cocking Foot Control - Standard when Auto-Cocking Drops are ordered
- Single Foot Pedal Elevation 120V (Plastic)
- 23" width seamless cushions with multi-density foam
- Paper Roll Holder and Cutter
- Locking 6" Extending Ankle Rest
- Lifting Capacity: 400 lbs

### Pelvic Options

- ES95009 - Elevating Pelvic Section w/single side control

### Elevating Foot Control Options

- 95284 - Dual Foot Pedal Elevation (Plastic)

### Upgrade To Auto-Cocking

- ES9320 - Elevating Table w/Auto-Cocking  
The patented Accelerator III Auto-cocking drops include:  
Thoracic, Lumbar and Pelvic Drops. (Patent# 7,381,214B1)  
\*Note: Air Compressor Sold Separately

### Air Compressor

- E7069 - Super Silent 20 Amp Compression - 115V
- E7068 - Super Silent 20 Amp Compression - 230V

### Headpiece Options

- ES95002 - Tilt, Forward Drop, Adjustable Face Cushions
- ES95003 - Tilt, Elevation, Adjustable Face Cushions
- ES95004 - Tilt, Elevation, Forward Drop, Toggle Drop, Adjustable Face Cushions

### Chest Piece Options

- ES95005 - Breakaway, Elevating Chest, with dual side control

### Drop Options

- ES58353 - Thoracic Drop
- ES58354 - Lumbar Drop
- ES58355 - Pelvic Drop



See our website for complete table details



Conforms To CAN/CSA C22.2 No. 601.1  
IEC 60601-1-4 • IEC 60601-1 • UL STD 60601-1





“I appreciate the versatility of the ES2000. It descends low enough for side-posture manipulation and elevates high enough for other manual treatment and examination. Since clinic space is valuable, the versatility of the ES2000 allows me to perform a variety of manual treatments on one excellent piece of equipment.”

— **Bradford J. Cole, DC, MS, CSCS**  
*Cole Pain Therapy Group, Bartlett, TN*



Standard 6"  
Extending Ankle Rest



Patented Accelerator III  
Manual & Auto-Cocking  
Drops



ES95005 - Breakaway,  
Elevating Chest Option



ES95004  
Headpiece Option



ErgoStyle ES2000 GEN II  
part number ES9300  
shown with options

# ErgoStyle<sup>TM</sup> FX Counterbalanced Flexion Table

## World-class Patented DROPS, Upgradeable In The Field

### ErgoStyle FX - ES5820

Elevating Manual Flexion Table

Starting at \$6,158

### ErgoStyle FX - ES5822

Auto-Cocking w/3 Drops, Elevating Manual Flexion Table

Starting at \$9,468

ErgoStyle<sup>TM</sup> FX Counterbalanced Flexion addresses the core needs of chiropractors. The ErgoStyle<sup>TM</sup> FX table is a sophisticated table that offers full, ultra smooth, counterbalanced flexion at its lowest height (20") and quietly elevates to its full height of (31"). The new design offers infinitely variable electric spring tension distributed in such a way to allow for fingertip effort in moving the caudal section in any direction.



**NEW OPTION**  
Long Axis Distraction

ES58494 - Manual Flexion  
ES59496 - Auto Flexion  
Long Axis Distraction Option



ES95004  
Headpiece Option



ES58359 - Breakaway,  
Elevating Chest Option



ErgoStyle FX  
part number ES5820  
shown with options

### ErgoStyle FX Features

- Electric Elevation from 20" to 31" with single foot control and expanding covers over lifting mechanism
- Manual flexion/distraction with variable electric spring tension allows smooth, counterbalanced flexion, extension or circumduction with the touch of a finger (flexion 15°, extension 12°, side bending 40° total)
- Tilting headpiece with adjustable face cushions
- Modular design for future upgrades
- 23" width seamless cushions with multi-density foam
- Stationary Chest
- Single Elevation Foot Pedal (Plastic)
- Single Tension Foot Control - Manual Flexion Only (Plastic)
- Paper Roll Holder and Cutter
- Locking 8" Extending Ankle Rest with ankle cuffs
- Lifting Capacity: 400 lbs

See our website for complete table options available.



# ErgoStyle<sup>TM</sup>HYLO



## World-class Patented DROPS, Upgradeable In The Field

### ErgoStyle HYLO - EH9510

#### Standard HYLO Table

Starting at \$8,497

### ErgoStyle HYLO - EH9530

#### Auto-Cocking w/3 Drops Table

Starting at \$11,762

ErgoStyle<sup>TM</sup> HYLO is a revolutionary adjusting table with many standard features previously considered upgrades. Improved ergonomics include: bilateral control placement, powered front and pelvic sections, free falling Accelerator<sup>TM</sup> III Drops, gas spring assist on headpieces, chest and pelvic cushions.



ErgoStyle HYLO  
part number EH9510  
shown with options

### ErgoStyle HYLO Features

- Improved ergonomics that include bilateral control placement, handle design, free falling accelerating drops and gas spring assist on optional elevating head, chest, and pelvic sections.
- Constant flow components allow a constant tilt speed with Automatic Weight Compensation regardless of patient's weight and height
- Powered front and rear sections allow full power adjustment of table length
  - Powered chest extension: Extends 14.5"
  - Powered pelvic extension: Extends 14.5"
- Controlled safety system: 2 step tilt down control
- Patented Accelerator<sup>TM</sup> III Manual Cocking Drops (Patent# 6,679,905)
- Fixed table height: 23.5"
- Table length adjustment: 54" - 83"
- Tilting mechanism: Gear drive hydraulics with a "soft stop" in the down position
- Tilting headpiece with adjustable face cushions
- Modular design for future upgrades
- 23" width seamless cushions with multi-density foam
- Stationary Chest
- Fixed Pelvic Section
- Single side controls to power chest and pelvic extension
- Single Foot Tilt Controls/Hylo (Plastic)
- Paper Roll Holder and Cutter
- Footplate: Foldaway
- Fixed Ankle Rest
- Lifting Capacity: 400 lbs

See our website for complete table options available.

# DOC DECOMPRESSION

Has the ability to correctly position the patients in antalgic posture

## DOC DECOMPRESSION - E9011-B

Targeting of Axial Rotation and Lateral Flexion

90 Days of FREE clinical support from consultant, Dr. Bryan Hawley

Starting at \$12,999

Take your Chiropractic practice to the next level with Pivotal Health Solutions DOC Decompression Table. The unique design of the new DOC table utilizes a state-of-the-art digital command center for specific vertebral targeting, separate lumbar and cervical decompression programming and continuous readout and graphing of treatment protocols. In addition, the digital command center also controls table elevation and rear flexion, all at the push of a button.

## DOC Features

- Axial Rotation and Lateral Flexion can be done simply with the press of a button and manipulation of the guide bar at the rear of the pelvic section. The doctor or practitioner can move the flexing section to its proper position such as lateral flexion or axial rotation. This is unlike other units that have a direct and straight pull only. In addition, this system comes with P-A flexion capabilities.
- Digital Command Center which includes:
  - Specific vertebral targeting
  - Separate lumbar and cervical decompression programming
  - Continuous readout and graphing of treatment parameters
  - Interactive color display
- Lifting Capacity: 350 lbs
- Pre-programmed lumbar and cervical decompression protocols
- Cervical capture headpiece
- Customizable treatment protocols
- Traction belting system can enhance restraint for optimal decompression.
- Reference library - including patient setup and training materials
- User-defined treatment hold and relax times
- Speed Control
- Expanded customizable protocols
- Real-time digital treatment tracking
- Table elevation
- Comfort cushions

“My practice, treats exclusively acute and chronic disc pathology. We consider ourselves to be the best in the field at what we do, and our patients count on us to deliver results. That is why after looking at all other decompression tables we have chosen the DOC decompression table by Pivotal Health. Not only can we treat our patients supine or prone but also from various angles such as flexion, extension, lateral side bending, and even axial rotation!”

— Dr. Michael L. Hall, D.C.  
Triangle Disc Care, Raleigh, NC



Cervical  
Capture Headpiece

E7100  
DOC Strapping System

See our website for complete  
table options available.



# HTG DROP Bench

Built from solid hardwood and quality materials



## HTG - HTG100 HTG Three Drop Bench

Starting at \$1,425

## HTG - HTG200 HTG Roller Massage Bench

Starting at \$950

Your new HTG Series has been designed and engineered to comply with the rigid standards of practitioners. The HTG 100 bench and HTG 200 Roller Massage Bench come with many standard features. Your new bench is designed to provide you with years of service. With proper maintenance your bench should remain virtually trouble free.



HTG -100  
shown with options



HTG -200  
shown with options

### HTG -100 Features

- Tilt headpiece - allows 30° of forward flexion and 30° of extension
- Five heights: 18", 20", 22", 24", 26"
- Cervical drop
- Pelvic drop
- Thoracic drop
- Drops are a full 0.75"
- Beautiful solid oak varnished legs
- Thoracic and lumbo-pelvic drop handles are on both sides of the table with the tension control on the left. The cervical drop handle and tension control is found on the right side of the headpiece.
- 2.5" firm foam on the main table body cushions, providing a very supportive adjusting surface. Headpiece cushions use 2.5" medium foam for patient comfort.
- Paper holder and cutter - accommodates standard 9" by 125' rolls
- Arm rests - 3.75" x 12". Arm rest foam is 1.5" firm foam
- Capacity: 450 lbs

### HTG -200 Features

- 8 rollers on a 24" track
- Spring loaded carriage
- Built in foam pad
- Extra wide 24" table x 70" long
- Height available from 17" - 30"
- Face hole cutout
- 30 Minute timer
- Beautiful solid oak varnished legs
- Capacity: 400 lbs

See our website for complete table options available.

# ErgoStyle<sup>TM</sup> Benches

Quality, Ergonomics, Performance. It's all in the family



EB9020  
ErgoBench Fixed Top  
Soft Foam



EB9000  
ErgoBench Tilt Headpiece  
Soft Foam

## ErgoBench - EB9020/EB9030

Fixed Top - Soft Foam/Firm Foam

Starting at \$668

## ErgoBench - EB9000/EB9010

Tilt Headpiece - Soft Foam/Firm Foam

Starting at \$778

The ErgoBench brings uncompromised stability and style to any office at an affordable price. It has an 18" to 24" height adjustment and features a unitized steel frame with a bottleneck shape, seamless cushions, multi-fit arm rests, and a paper roll holder and cutter. It is also available with a comfortable and easily adjustable tilt headpiece. Ships with unassembled legs. The ErgoBenches can be shipped completely assembled for an additional charge.

## ErgoBench Features

- Unitized steel frame
- Seamless cushions
- Multi-fit arm rests
- Paper roll holder and cutter
- Bottleneck shape
- Available in soft or firm foam
- Height Range 18" - 24"
- Capacity: 350 lbs



EB9040  
ErgoBasic Tilt Headpiece  
Soft Foam



EB9040  
ErgoBasic Soft Foam  
with optional ES95004  
Headpiece

## ErgoBasic - EB9040/EB9050

Tilt Headpiece - Soft Foam/Firm Foam

Starting at \$2,262

A critical mix of adjusting features that is praised by many practitioners. A tilt head and accelerator pelvic drop are standard features. Optional extras include an upgradeable headpiece, allowing operation of tilt, adjustable face opening, forward drop, elevation and toggle drop. See complete headpiece options below. Gas spring powered pelvic elevation is also an option, allowing you to effortlessly put patients into passive flexion. ErgoBasic is sleek, elegantly simple, yet affordable.

## ErgoBasic Features

- Tilt headpiece
- Accelerator<sup>TM</sup> III pelvic drop
- Paper roll holder and cutter
- Unitized steel frame
- Seamless cushions
- Available in soft or firm foam
- Shipped fully assembled
- Height Range 18" - 24"
- Capacity: 350 lbs



## ErgoBasic Options

- ES95001 - Tilt headpiece with Adjustable Face Cushions
- ES95002 - Tilt headpiece w/Forward Drop and Adjustable Face Cushions
- ES95003 - Tilt headpiece w/Elevation and Adjustable Face Cushions
- ES95004 - Tilt headpiece w/Elevation, Forward Drop, Toggle Drop, and Adjustable Face Cushions
- EB89295 - Elevating Pelvic Section



# BenchMark Therapy Tables

The Benchmark Series of exam/therapy tables delivers uncompromising durability, versatile features and affordable prices to any office looking for an entry level solution

## Single Section - E9139

Therapy/Exam Table (28"x78"x30")

Starting at \$715

## Double Section - E9143

Therapy/Exam Table (28"x80"x30")

Starting at \$818

## Triple Section - E9142

Therapy/Exam Table (28"x80"x30")

Starting at \$921

### BenchMark Features

- Heavy duty steel frame
- Durable powder coat finish (black only)
- Seamless cushions
- 2" high density foam thickness
- Paper roll holder and cutter
- 28" wide cushions
- High impact levelers
- Face slot
- 30" fixed table height
- 18" storage shelf
- Capacity: 350 lbs

\*Note - Ships with unassembled legs. The Benchmark benches can be shipped completely assembled for an additional charge.

E9139



E9143



E9142



# Clinic Rolling Stools

## Rolling Stool 19332 Features

- Working weight: 300 lbs
- Seat padding: 4" deluxe cushion
- Height Adjustment: 19"- 24"
- Five star base with heavy-duty casters.
- Pneumatic controls

\$65



19332  
Rolling Stool

## Rolling Stool 19310 Features

- Working weight: 300 lbs
- Seat padding: 4" deluxe cushion
- Height Adjustment: 19"- 24"
- Back Height Adjustment: 30"- 35"
- Removable back support
- Five star base with heavy-duty casters.
- Pneumatic controls

\$75



19310  
Rolling Stool with Back

# Intersegmental Roller Tables

Pretreatment Therapy and massage tables



ATT-300  
Standard 22"  
wide top



ATT-300  
Optional 25"  
wide top shown  
with cervical and  
knee bolsters

## ATT-300

Roller Massage Table (72"L x 22"W x 25"H)

Starting at \$2,466

The ATT-300 Intersegmental Roller Massage Table was specially designed for pretreatment therapy and massage. Under its thick, padded lounge top, three large rotating rollers travel back and forth, kneading the muscles in the neck, shoulders and lower back.

### ATT-300 Features

- Clockwise and counter-clockwise rotation
- Travel indicator
- Automatic lowering device
- Adjustable massage pressure
- Vibration
- LED roller height indicator
- 30-minute timer
- Cervical and knee bolsters
- Roller height adjustment 4"
- Treatment area 25"
- Capacity: 350 lbs



ErgoWave shown  
with optional knee  
bolster

## ErgoWave - EW9080

Roller Massage Table (78"L x 25"W x 24"H)

Starting at \$2,497

Unique Wave Action – The ErgoWave roller pattern travels up the spine like a wave. All areas of the patient's spine are covered with three independent rollers repeating every 13 seconds.

### ErgoWave Features

- Triple Rollers mounted on three axles to isolate lumbar, thoracic and cervical regions
- Quick Height Adjustments– A gas spring makes roller height adjustments quick and easy
- Ultra Quiet Operation - High quality motor for nearly silent operation without annoying vibrations
- Built-in storage shelf
- A welded steel frame ensures stability
- Roller height adjustment 3"
- Treatment area 25"
- 30-minute timer
- Capacity: 350 lbs



Conforms To CAN/CSA C22.2 No. 601.1  
IEC 60601-1-4 • IEC 60601-1 • UL STD 60601-1

# Portable Tables

Built from quality materials to withstand rigorous conditions

## Chiron Pro Stationary - I-ChironPro-PU 2 Drop Fixed Table Height 19"

Starting at \$822

### Chiron Pro Stationary Features

- 1.5" Chiropractic Foam System
- Five Position Headpiece
- Thoracic and Lumbar Drops
- Hand Rests
- Paper Holder
- Carry Case
- Table Dimensions: 66" L x 21" W x 19" H
- Capacity: 400 lbs
- Table Weight: 43 lbs

The Chiron Pro Stationary offers a five position headpiece with Lumbar, and Thoracic Drops, carry case and a fixed height of 19". This table is tough enough to withstand the requirements of the most demanding use. \*Imported



## Chiron Pro Adjustable - CHID 2 Drop Adjustable Height Table

Starting at \$1,234

### Chiron Pro Features

- 1.5" Chiropractic Foam System
- Five Position Headpiece
- Thoracic and Lumbar Drops
- Hand Rests
- Paper Holder
- Carry Case
- Table Dimensions: 66" L x 21" W
- Adjustable Legs 19.25" to 25.25" (with height adjustments every 2")
- Capacity: 400 lbs
- Table Weight: 43 lbs

We took the dependable Chiron Table one step further by adding smoothly operating lumbar and thoracic drops and standard hand rests. By truly addressing the need for a portable table that is uncomplicated and durable, yet advanced enough to accommodate a variety of adjustment techniques, the Chiron Pro has quickly become the table of choice for today's chiropractor on the go. MADE IN THE USA



## Omni Sideline - OM-3007 800 lb Capacity - Portable Treatment Table

Starting at \$599

### Omni Sideline Features

- 3 pocket carry case
- Length: 73"
- Height Range: 25"-35"
- Width: 28", 30"
- 2" Firm foam system
- Removable field feet
- Dual headrest ports
- Inverted Truss™ system
- Capacity: 800 lbs
- Table Weight: 37 lbs

For the mobile trainer/Chiropractor, the Omni Sideline table is lightweight and tough. This portable treatment table features solid wood construction, 2 inch high-density foam top and a folding design to allow for easy storage, set-up & portability







"Apollo is by FAR the best designed, most efficient and easiest to use laser I have ever seen. It is durable, reliable and therapeutically superior. My patients love it, my practice and income have grown from it, and I can't imagine how I could practice without it. If you compare the features, design, trustworthiness, and VALUE of low level lasers available today, there is no better choice than the Apollo. "

— **Dr. Mark de Dubovay**  
*Director, Advanced Wellness Center  
 Long Beach, CA*

## apollo Laser Therapy

### Introducing the Newest Line of Apollo Lasers

**NEW**  
**Handheld**  
**2,000 mW LASER**  
**- HH2000 -**

**Apollo 2,000 mW Cold Laser**  
 Low Level Laser therapy just became more flexible. Experience the freedom of this incredible new device.

Starting at \$3,995



Small, lightweight unit fits in the palm of your hand. Simple on and off button, no programming necessary. The new HH2000 handheld device from Apollo uses low level laser technology to speed the recovery process. The HH2000 laser is a non-invasive, safe way to manage pain—without negative side effects. Includes: unit docking station (shown), carry case, 2 safety goggles, battery pack (shown), belt holster and everything you and your patient need for a successful treatment.

The Apollo laser efficiently harnesses advanced Class 4 laser technology. Using Apollo's perfect marriage of power and wavelength, you are assured fast, safe, comfortable treatments, exceptional clinical results and an impressive bottom line. And unlike hot Class 4 lasers, the Apollo is cool, safe, and effective.



**90 Days FREE Clinical Support**  
 from renowned speaker, educator, and author - **Dr. Curtis Turchin**

Pivotal Health has strategically brought together some of the industry's leading consultants in their respective fields in order to support your success. With the purchase of any Apollo Laser System you will receive 90 Days of FREE Support from Dr. Curtis Turchin



## Apollo Desktop Laser Unit

### 2 Channel Desktop with Choice of 500mW, or 3,4 or 5,000mW Probe

Starting at \$5,897

- Heavy-duty construction for office and clinic environments
- Probe holders on both sides of the unit that can accommodate two probes and both right and left handed practitioners
- Elevated face for easy viewing and control
- 2 timers allow you to treat 2 patients at the same time

## Apollo Portable Laser Unit

### Portable Control Unit with Choice of 500mW or 3,000mW Probe

Starting at \$6,051

- Weighing only 2.3 pounds, the Apollo laser is the lightest, most powerful portable laser you can buy
- Custom Carrying Case allows for easy travel from office to office and out in the field
- User replaceable battery
- Heavy-duty construction for office and clinic environments

See our website [www.ApolloPT.com](http://www.ApolloPT.com) for details

## Desktop & Portable Laser Unit Features

- LCD display that provides probe status and treatment times
- User controlled treatment times selectable between 10 seconds and 2 minutes
- Built in safety and fault detection software
- Built in power test for assessing probe output
- Durable aerospace quality aluminum assembly
- ISO 13485 quality assured manufacturing

## Desktop & Portable Laser Units Include

- APOLLO Control Unit, Desktop or Portable
- Choice of (1) laser probe: 500mW, 3, 4 or 5,000mW
- Power Adapter
- (1) Probe Cable
- (2) pair of Safety Goggles
- Laser Carrying Case
- Book: Light & Laser Therapy: Clinical Procedures
- Patient education brochures
- Laser tutorial DVD with systems



## 3,4 & 5,000mW Point Probes

### 810nm Infrared Cluster probes

3,000 mW (4x750mW)  
4,000 mW (4x1000mW)  
5,000 mW (4x1250mW)  
Apollo uses premium medical grade Class 4 laser diodes that provide extra long life and reliability. The probes are made from aerospace grade aluminum which means cool durable performance.



## 500mW Point Probe

### Hand, Small Area Laser

The 500mW point probe was designed specifically to reach those difficult to reach areas or areas with less tissue to penetrate such as hands, feet, trigger points and auricular areas.

Optional fine tip shown.



# Whole Body Vibration Therapy

## Vertical Movement and Whole Body Vibration




V1000 shown

Whether you want to reduce weight, improve flexibility or even improve your golf swing, the new A3 Vibration will help you reach your goals. Vibration training is low impact and gentle, yet unbelievably effective. Anyone, regardless of their age or fitness levels, can now benefit from exercising without straining their joints or other tissues.

A simple touch of a button and the A3 Vibration unit turns into a pleasant vibrating platform. This then transfers the vibrations to muscles in the body causing them to react gently yet powerfully in a reflex action.

By standing on the unit you will notice how your body automatically adjusts to the vibrations. This mechanical stimulus produces a stretch reflex which (vibratory reflex), depending on the selected frequency, results in muscles gently contracting 20 to 60 times per second. Utilizing the body's innate reflexive response to disruptions in stability, Whole Body Vibration stimulates and enhances muscle strength and performance.

Different parts of the body are put in line with the vibrations. The user then assumes various exercise positions and varying intensities, ranging from squats, lunges and push-ups to abdominal crunches, tricep dips or bicep curls. Poses are typically held for just 30 seconds to 1 minute at a time, during that timeframe, multiple muscle groups are activated simultaneously to stretch, tone and strengthen the target tissues.



## Welcome Clinical Consultant Christian Reichardt, DC



Christian H. Reichardt has been a practicing doctor of chiropractic since 1983. He was a disability evaluator for the State of California, but then turned his attention from work-related injuries to sports injuries due to his own involvement in sports. He became one of the first Certified Chiropractic Sports Physicians in the USA in 1986, as well as a Diplomat of the American Academy of Pain Management. His practice in Santa Monica, California draws a large clientele of Hollywood celebrities, as well as professional and amateur athletes. He has worked with a number of US Olympic Team members and Golf professionals.



Scan to learn more





Dimensions:  
(L) 36.7in.  
(W) 29.5in.  
(H) 56.1in.

**27**  
Exercise  
Postures

\$1,999

The **V1000** is an easy to use vibrating platform that is equipped with a user-friendly virtual coach console. Mainly for lower body training but can help upper body and core training. The training will improve muscle endurance, explosive strength and bone density while helping to avoid exercise injuries and muscle atrophy. It also can be used as an aid for warm up or cool down prior or after exercising.

### V1000 Home & Small Clinic:

Features a small "footprint" to make it easy to fit in your space and programs that suit the majority of patients. Ideal for light duty clinical environments and home use.

- Medium-intensity Vibration Training
- 4 Training Programs : Total Body, Upper body, Lower body, Abdominal
- 1 Manual mode : user can adjust time & frequency
- Workout Display: Posture Exercise Program, Exercise Time, Frequency, Exercise Group
- Frequency: 20 Hz ~ 50 Hz
- Amplitude: 1.5 mm
- Force: 800 N
- Maximum Load: 308 lbs
- Dimension: L 36.7" x W 29.5" x H 56.1"



Dimensions:  
(L) 39in.  
(W) 38.3in.  
(H) 63.8in.

**64**  
Exercise  
Postures


\$3,999

The **V3000** is an exceptional vibration platform for its multi-purpose functions. Besides normal vibration benefits, the V3000 has a built-in professional sculpture fitness program. With this program, users not only train their body but also sculpt their shape all on this single machine. This platform supports user flexibility and stability.

### V3000 Gym & Large Clinic:

Features a larger platform and more programs for the largest range of exercises and patient types.

- Medium-intensity Vibration Training
- 8 preset program models, 10 sets of editing exercise models & 1 manual model
- Training Programs:  
WHOLE BODY - Strength & Sculpting  
UPPER BODY - Strength & Sculpting  
LOWER BODY - Strength & Sculpting  
CORE TRAINING - Strength & Sculpting
- High and low force of vibration.
- Frequency: 25 Hz ~ 50 Hz
- Amplitude (MAX: HIGH: 2.5mm; LOW: 1.2mm)
- Force: HIGH: 2.5KN; LOW: 1.8KN
- Maximum Load 350 lbs
- Dimension: L 39" x W 38.3" x H 63.8"
- Heavy Duty, Commercial Grade 1 HP Motor



### All A3 Vibration Units Feature:

- A backlit LCM console display with clearly visible workout status
- User-friendly virtual coach console
- Non-slip handlebar allows easy use & training from various angles
- Unique dual tensile-straps offer user enhanced training for upper arms
- Wide vibrating platform with a stabilizing and anti-seismic rubber mat
- All units come with rubber mat, and exercise/stretching straps
- Exercise protocols and stretching posters
- Voltage: 120 VAC / 60 Hz
- Motor: DC Motor 1.0 HP

# Teeter Inversion Therapy

## Certified Quality Inversion Equipment CE UL



Teeter EP-550  
shown

### Teeter EP-550 - EP1001

#### Base Model Manual Home-Use Inversion Table

Starting at \$329

You've seen us on TV, you've heard of your favorite athlete using a Teeter... your neighbor even has one. Now it's your turn to try Teeter Hang Ups, the inversion market pioneer since 1981. With the passion of someone who has lived through back pain and found relief, Roger Teeter has built Teeter Hang Ups into a brand recognized worldwide for superior quality, function and innovation.

### Teeter EP-550 Features

- Ergonomic Flex-Technology™
- Optimum Traction Surface
- No Tear, No Stain Construction
- Ankle Comfort Dial™
- Precision Balancing
- Pre-set Rotation Control
- Customize With Accessories
- Stable Base
- Unique and Patented Security Features
- Smooth, Quiet Rotation
- Support Hand Grips
- Quick Storage



Teeter EP-950  
shown

### Teeter EP-950 - EP1019

#### High End Manual Home-Use Inversion Table

Starting at \$399

Teeter has exclusively designed and manufactured inversion equipment that adheres to and proudly displays the UL and CE quality listing and as a result has built a reputation for incomparable quality standards, testing and marketing efforts. Teeter inversion products give healthcare providers a solution for all of their back pain patients. Teeter's complete product offering allows clinicians to treat patients in the office, clinic or training room and then, by use of a Teeter home inversion table, manage the patients back pain outside of the office.

### Teeter EP-950 Features

- Ratchet Ankle Lock System
- Ergo-Embrace Ankle System™
- EZ Stretch Traction Handles™
- Over EZ Handles™
- Pre-set Rotation Control
- Robust Medical-Grade Frame
- Ergonomic Flex-Technology™
- Optimum Traction Surface
- No Tear, No Stain Construction
- Ankle Comfort Dial™
- Precision Balancing



### Teeter DFM - M1-1020

#### Medical Prone/Supine Power Inversion Table

Starting at \$3,545



Teeter clinical products, the DFM and DEX inversion units are used by chiropractors, physical therapists and athletic trainers at the professional and collegiate level. Experience the versatility of the DFM by maximizing function in a small space. The DFM is an effective therapeutic tool for any clinic. Patients may rest prone or supine while inverted, allowing for treatments at varied angles. Adjustable features ensure user comfort while the easily accessible push-button control allows for supervised or unassisted inversion sessions.

See our website for complete listing of Teeter Units available.



# FMST Tools


## Fluid Motion Soft Tissue Tools




### Welcome Clinical Consultant Jerad Bludorn, DC

Dr. Jerad Bludorn, DC is a Chiropractic Physician and Acupuncturist with an extensive background in soft tissue therapy, functional movement and stability. He has trained extensively with Motion Palpation Institute as well as many other treatment techniques and world renowned lecturers and researchers. He has over 10 years experience in private practice working with Orthopedic Surgeons and Medical Doctors to improve outcomes for his patients, and developed FMST Tools to assist him in improving his outcomes. Dr Bludorn has successfully implemented FMST Tools into many different types of practices including Physical Therapy, Chiropractic, Massage Therapy, Occupational Therapy and Athletic Training.

See our website for complete details



F1 shown



F2 shown



### FMST Tools

#### Fluid Motion Soft Tissue Tools - (set of 2)

Starting at \$699

FMST Tools are the soft tissue mobilization tools that can take myo-fascial manipulation and mobilization to the next level in your practice. The design of the (2) tools available enables you, the practitioner, to work over all surfaces of the body quickly and easily. We have been able to reduce scar tissue problems from old injuries/surgeries and achieve greater results with pain control, joint and muscle flexibility that many thought was unattainable.

- 2 tool designs to treat all areas of the body.
- Unique metal injected plastic material for superior weight and durability over other tools.
- Custom molded non-slip thumb grip to create maximum angular pressure without over use by the practitioner.
- Soft but aggressive treatment edge to work fascial adhesions, post surgical scarring and muscle injuries.
- Ambidextrous Design: No more switching through multiple tools to find the right edge or angle. Both tools are designed to be used on both sides and on every surface for increased compatibility of treatment on multiple body areas.



# Pressure Positive Body Tools

## Self Care Deep Muscle & Injury Prevention Body Tools

### The Original **Palmassager®**

One of the most efficient massage therapy aids available anywhere. Hold it firmly in the palm of the hand and use one, two or three of its round knobs to massage sore muscles.



~~\$9.95~~ **\$7.95**

### The Original **Index Knobber II®**

A simple device designed to let you apply deep pressure to muscles within easy reach of your hand. The ideal instrument for clinic or home use allowing for firm, sustained pressure on targeted muscles.



(Wood option)

Index Knobber ~~\$9.95~~ **\$7.95**

Wood Index Knobber (w/ ball) ~~\$9.95~~ **\$7.95**

### The Original **Knobble II®**

The new design of the all wood favorite. Made of durable polymer with a band of non-latex Santoprene, The Knobble II gives your hand more power when you need it.



Knobble II . . . . . ~~\$9.95~~ **\$7.95**

Wood Knobble . . . ~~\$17.95~~ **\$15.95**

### The Original **Vibrassager®**

This innovative rolling massager is equipped with two knobbed wheels that actually feel like they vibrate as they roll across the back giving a delightful, relaxing massage.



~~\$17.95~~ **\$15.95**

### The Original **Orbit Massager®**

The Orbit Massager® allows you to provide a stress releasing massage to the back, neck, shoulders and muscles throughout the body.



~~\$24.95~~ **\$18.95**

### The Original **Jacknobber II®**

An elegantly designed hand-held massage tool that can be used to provide a more pleasant and relaxing massage.



~~\$9.95~~ **\$7.95**

### **Tiger Tail®**



Make your muscles happy! The Tiger Tail Rolling Muscle Massager, 18" or 22" is a dense foam covered, muscle-friendly self care massage tool that helps you apply rolling and compression massage techniques.

18" - TT18 ~~\$29.99~~ **\$24.99**

22" - TT22 ~~\$34.99~~ **\$29.99**

### The Original **Backnobber II®**

The Original Backnobber® II is constructed of durable, fiberglass reinforced Nylon and carries a full satisfaction guarantee.

The Original Backnobber® II is used by hooking one of its curves over a shoulder or under an arm and using it's leverage to apply firm, deep pressure to knots in the muscles of the neck, shoulder, upper and lower back and to nearly any other location where muscle soreness may occur.

~~\$41.95~~ **\$31.95**



**FREE** 35page  
fully illustrated User Guide.



ReSell Program Available - Call 800-627-2387 for details

### Large Floor Display

This flexible, 2-sided rolling display unit is built from solid wood and will hold a maximum of 81 pieces. You can mix and match Pressure Positive products on this unit to best fit your needs.

**FREE Display**

~~\$1,649.15~~ **\$1,300**  
Free display with purchase  
of 81 tool set

Products	Qty
Knobble II	17
Backnobber II w/Manual	20
Vibrassager	4
Indexknobber II	12
Palmassager	8
Orbit Massager	3
Jacknobber II	12
18" Tiger Tail	3
22" Tiger Tail	2

### Counter Display

This flexible counter display is built from solid wood and will hold a maximum of 36 pieces. You can mix and match Pressure Positive products on this unit to best fit your needs.



Products	Qty
Knobble II	4
Backnobber II	8
Vibrassager	4
Indexknobber II	4
Palmassager	8
Orbit Massager	4
Jacknobber II	4

~~\$926.24~~ **\$775.24**  
display sold separately



# RelaxRight Pillows

Dr. Designed Orthopedic Memory Foam Pillows

Made with  
**Memory Foam**



## Adult Orthopedic Memory Foam Pillow

**AWARDED 21 Patent Claims!**

Adult Standard 90103 - \$38.50 / MSRP \$59.99

Unlike any other pillow on the market, the adult Relax Right® pillow was designed by a doctor and was awarded 21 patent claims. Relax Right® combines two types of high quality memory foam and a hollow core to deliver the ultimate pillow for both back and side sleeping. 100% Hypoallergenic. Antimicrobial velour cover. 23" x 13" x 4.5"



## Adult Contour Memory Foam Pillow

**Luxurious comfort and support**

Adult Contour 90102 - \$32.50 / MSRP \$49.99

The contour-shaped pillow is made of temperature-sensitive memory foam that molds to your body to provide luxurious comfort and support. The antimicrobial cover is made of plush cotton velour and is resistant to allergens, dust mites, mold and mildew. 24" x 12" x 4 1/2"



## Adult Ventilated Memory Foam Pillow

**Cool Circulating Core for exceptional Comfort**

Adult Ventilated 90101 - \$32.50 / MSRP \$49.99

The "soft feel" traditional-shape memory foam pillow is made of temperature-sensitive memory foam that molds to your body to provide luxurious comfort and support. The cool circulating core helps to eliminate overheating. The antimicrobial cover is made of plush cotton velour and is resistant to allergens, dust mites, mold and mildew. 24" x 16" x 5"



Big Kid shown

Toddler shown

## Childrens Orthopedic Memory Foam Pillows

**Recommended for ages 2-5 years and 5+**

Toddler Pillow 81003 - \$17.99 / MSRP \$29.99

Big Kid Pillow 81003 - \$25.25 / MSRP \$39.99

Did you know that ordinary pillows push the spine out of alignment? That's not good for your child's growing bones! Our RelaxRight® Toddler and Big Kid's Pillows position the neck and spine correctly, promoting healthier posture and better sleep. Our children's pillows are made with memory foam and feature an antimicrobial velour cover which adds to the great performance of our pillows! 15" x 10" x 3".

See our website for complete listing of RelaxRight Pillows available.

# Providing Solutions For Five Markets

Chiropractic | Massage | Therapy | Medical | Athletics



**PHSTherapeutics**  
by Pivotal Health Solutions

## Manufacturing the premier Functional Rehabilitation table

Pivotal Health Solutions, a leading manufacturer of premium quality health and wellness equipment for over 20 years, proudly unveils their newest brand, "PHS Therapeutics". Pivotal Health Solutions has earned an excellent reputation in the markets they serve and believe this is why their products are used by so many of the world's most prestigious schools, professionals, teams and therapists. In addition to their dedication to quality, many of their customers value the unique features available only from Pivotal Health Solutions. *"We continue to strengthen our company through organic growth, acquisitions and strategic relationships. The addition of the PHS Therapeutics brand of leading line Therapy Products broadens our product offering and supports our initiative to more effectively support Therapy Professionals," said Jim Engels, CEO of Pivotal Health Solutions.*

PHS Therapeutics offers HyLo Electric Treatment Tables, high quality solid Oak Wood Treatment Tables, Work Assist Stations, Vibration Therapy, Inversion Therapy, Cold Laser Therapy by Apollo, the patent pending (LAST) Leg and Shoulder Therapy Table, Mat Tables, Mobile carts, Storage Cabinets, Medical Cabinets, Weight Carts, (FMST) Fluid Motion Soft Tissue Tools, along with a full line of Portable Treatment Tables.



**CUSTOM  
CRAFTWORKS**

A Brand Of Pivotal Health Solutions

## Providing Best-In-Class Massage Tables for over 28 years

Custom Craftworks launched in 1986, carving out a niche for high-end, quality manufactured, customized massage tables that meet the ergonomic needs of therapists and the comfort requirements of clients.

Since then, the company has grown to become one of the largest manufacturers of massage tables, chairs, equipment and accessories in the United States. Acquired by Pivotal Health Solutions in 2009, today Custom Craftworks continues to build thoughtfully constructed, long-lasting products for students, therapists and educators in the massage, therapy and holistic health fields.

The company's accessories line also meets these high standards for quality and innovation. Custom Craftworks is known for developing a number of specialized accessories that are now staples in the field. Our product line includes the popular Sani-cover™ for face rests, to massage therapy self-care tools, and Omni™ cervical relief pillow. These are products that both clients can enjoy and massage therapists can utilize in their practice.





# Warranty / Upholstery

## Stock Vinyl Upholstery Colors UPHOLSTERY CARE AND CLEANING GUIDE

For light soiling, use a solution of 10% household liquid dish soap with warm water, applied with soft damp cloth. Rinse with clean water and dry. For heavy soiling, dampen a soft white cloth in a one to one (1:1) solution of Fantastik® and water OR Formula 409® and water. Rub gently and rinse with a water dampened cloth. For more difficult stains, dampen a soft white cloth with a solution of household bleach (10% bleach / 90% water). Rub gently and rinse with a water dampened cloth to remove bleach concentration. **DO NOT USE ALCOHOL BASED CLEANING AGENTS! IMPROPER CLEANING WILL VOID THE WARRANTY.** This information is not a guarantee and does not relieve the user from the responsibility of the proper and safe use of the product and all cleaning agents. The use of certain agents can be harmful to the surface appearance and lifespan of vinyl. Pivotal Health Solutions assumes no responsibility resulting from the use of such agents to the vinyl.



## Warranty Overview [see website for complete warranty details](#)

Pivotal Health Solutions warrants to the original retail purchaser that it will repair or replace components of the domestic and international products manufactured by PHS (except for components not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. Pivotal Health Solutions obligation under this warranty is limited to the repair or replacement, at Pivotal Health Solutions option, of the applicable components. This limited warranty shall only apply to defects that are reported to Pivotal Health Solutions within the applicable warranty period and which, upon examination by Pivotal Health Solutions, prove to be defective. This warranty extends only to the first retail purchaser of a product, and is not transferable or assignable.

## CHIROPRACTIC / THERAPY PRODUCTS WARRANTY

Products	Frame, Materials & Workmanship	Electrical & Moving Parts	Foam	Upholstery	Rollers
Whole Body Vibration	5 year parts and labor	2 year	N/A	N/A	N/A
ATT300, ErgoWave	1 year parts only, labor not included	1 year	1 year*	1 year*	Lifetime to original purchaser
ErgoBench, ErgoBasic	1 year parts only, labor not included	1 year	1 year*	1 year*	N/A
BenchMark Series	1 year parts only, labor not included	1 year	1 year*	1 year*	N/A
DOC Decompression	1 year parts and labor	1 year	1 year*	1 year*	N/A
TradeMark & TradeFlex	1 year parts and labor	1 year	1 year*	1 year*	N/A
ErgoStyle FX, ES2000, Hylo	1 year parts and labor	1 year	1 year*	1 year*	N/A
Portable Chiropractic Table	1 year parts and labor	1 year	1 year*	1 year*	N/A
Teeter DFM, DEX	2 year parts only, labor not included	2 year	2 year*	2 year*	N/A
Teeter EP550, EP950	5 year parts only, labor not included	5 year	N/A	N/A	N/A
Apollo Lasers	2 year unit parts, probes, and cables 1 year labor	Out of box warranty on battery	N/A	N/A	N/A
RelaxRight Pillows	N/A	N/A	N/A	N/A	N/A
FMST Tools	N/A	N/A	N/A	N/A	N/A

\*normal wear/tear excluded.

# PROUDLY KEEPING JOBS IN AMERICA

PHS Chiropractic, by Pivotal Health Solutions, was formed with the mission to be the leading manufacturer of quality built, ergonomically sound chiropractic adjusting tables. The growth of our trusted brand began with the launch of Eurotech Tables in 2005 and the acquisition of the Chattanooga ErgoStyle table line in 2008. Over the course of nine short years, Pivotal Health Solutions Chiropractic division has blossomed into one of the most comprehensive, quality built and ergonomically sound chiropractic adjusting table lines in the market today.



All PHS Chiropractic Tables are proudly engineered, custom made and built in America from US and globally sourced components.



page 19



page 16-17



page 18



page 8



page 10-11



page 12