

**Fluid UBE E920**

ITEM: E920



THE NEW, UNIQUE AND TECHNOLOGICALLY ADVANCED E-920 UBE has multiple adjustments to perfectly accommodate all levels of able bodied and those users with disabilities. The E-920 UBE Cycle is our top of the range medical grade rehab machine. The Upper Body Ergometer gives users the option of adjusting the crank arm length. This feature is especially useful for mechanical/ rehab/ and elite athlete use. The E-920 is perfect for the professional facility where maximum versatility of equipment is essential and is ideal for rehabilitation.

Key features include:

- 20 level patented Variable Fluid Resistance
- Variable length crank arms - ideal for mechanical/rehab use
- Instantly adaptable between 3 functions - seated, standing or wheelchair access
- Operates in forward & reverse directions
- Angle adjustable, self leveling hand grips
- Main arm adjusts from horizontal to vertical with internal gas-assist mechanism
- Interactive Performance Monitor
- 4 position, horizontal seat adjustment
- Contoured "soft touch" seat with lumbar support
- Foot rest for added stability
- Heart rate reception built in
- Compact footprint and self-generating
- ADA compliant
- Light Commercial Warranty: 10 years on the frame; 3 years Tank & seals; 2 years non-wearing components; All other 1 year